



**FLOW PLUS**  
DIRECTIONS FOR USE:

**A flowable light cure hybrid resin restorative suitable for Class III, Class IV, and Class V restorations**

FLOW PLUS is a flowable light cure hybrid resin restorative suitable for Class III, Class IV, and Class V restorations. Additional functions include: repair of enamel defects, repair of porcelain restorations, minor occlusal build-ups in non-stress bearing areas, pit and fissure sealant, cement for ceramic/composite veneers, incisal abrasions, and core build-ups. This unique composite which incorporates BIS-GMA chemistry along with a proprietary glass filler (approximately 60%) yields a flowable non-slumping restorative material.

FLOW PLUS combination of flowability and direct application system simplifies material placement and minimizes finishing to result in a consistently superior restoration.

**Directions for Use**

1. Isolate and dry preparation.
2. Etch enamel and/or dentin for 15 seconds with Medicept's Solo Etch (37.5% phosphoric acid gel\*). Rinse thoroughly and lightly air dry.
3. Apply a bonding system such as Bond Plus NT.
4. Select the appropriate shade. Dispense a layer of Flow Plus directly onto the preparation along the axial and gingival margins. Use a brushing action to allow the material to flow into the marginal interfaces and recesses of the preparation.
5. Incremental layers should be no more than 2 mm. Light cure each increment for a minimum of 40 seconds with a visible light curing unit.
6. Build contours and ridges to desired shape by flowing the composite into the preparation with a brush-like motion.
7. Finish and polish using conventional techniques.
8. For a glaze finish, apply one or two drops of Flow Plus, brush over the surface, and cure for 40 seconds.

• CAUTION: Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water, get medical attention for eyes. Do not take internally.

**Shading System**

Flow Plus is available in 16 Vita ccc shades and 5 special shades as follows:

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, Extra Light, G2 (gingival), Universal Opaque, White Opaque, Incisal,

**Notes:**

1. Replace syringe cap after each use to prevent the resin from hardening in the syringe.
2. Flow Plus dispensing tips are intended for single patient use only.

**Storage Conditions**

Store FLOW PLUS at room temperature. It has a shelf life of 2 years from date of manufacture.

(Note: VITA® is a registered trademark of VITA Zahnfabrik Company, Germany).

CAUTION: Uncured resin materials may cause contact dermatitis. Avoid prolonged exposure of uncured resins to skin and soft tissue.

**LIMITED WARRANTY**

**LIMITATION OF MEDICEPT'S LIABILITY**

Medicept's technical advice, whether verbal or in writing, is designed to assist dentists in using Medicept's product. The dentist assumes all risk and liability for damages arising out of the improper use of Medicept's product. In the event of a defect in material or workmanship, Medicept's liability is limited, at Medicept's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Medicept. In no event shall Medicept be liable for any indirect, incidental, or consequential damages.

**EXCEPT AS EXPRESSLY PRO-VIDED ABOVE, THERE ARE NO WARRANTIES, BY MEDICEPT, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.**

**THIS PRODUCT IS SPECIALLY FORMULATED FOR USE IN DENTISTRY**

**Medicept UK Ltd.**

2<sup>nd</sup> Floor, Hygeia House, 66 College Road, Harrow, Middlesex, HA1 1BE, United Kingdom.

Medicept Dental operates a policy of continuous monitoring and improvement of our products. If you have any comments about this product, please contact us at the above address stating the batch number shown on the packaging.



**FLOW PLUS**  
DIRECTIONS FOR USE:

**A flowable light cure hybrid resin restorative suitable for Class III, Class IV, and Class V restorations**

FLOW PLUS is a flowable light cure hybrid resin restorative suitable for Class III, Class IV, and Class V restorations. Additional functions include: repair of enamel defects, repair of porcelain restorations, minor occlusal build-ups in non-stress bearing areas, pit and fissure sealant, cement for ceramic/composite veneers, incisal abrasions, and core build-ups. This unique composite which incorporates BIS-GMA chemistry along with a proprietary glass filler (approximately 60%) yields a flowable non-slumping restorative material.

FLOW PLUS combination of flowability and direct application system simplifies material placement and minimizes finishing to result in a consistently superior restoration.

**Directions for Use**

1. Isolate and dry preparation.
2. Etch enamel and/or dentin for 15 seconds with Medicept's Solo Etch (37.5% phosphoric acid gel\*). Rinse thoroughly and lightly air dry.
3. Apply a bonding system such as Bond Plus NT.
4. Select the appropriate shade. Dispense a layer of Flow Plus directly onto the preparation along the axial and gingival margins. Use a brushing action to allow the material to flow into the marginal interfaces and recesses of the preparation.
5. Incremental layers should be no more than 2 mm. Light cure each increment for a minimum of 40 seconds with a visible light curing unit.
6. Build contours and ridges to desired shape by flowing the composite into the preparation with a brush-like motion.
7. Finish and polish using conventional techniques.
8. For a glaze finish, apply one or two drops of Flow Plus, brush over the surface, and cure for 40 seconds.

• CAUTION: Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water, get medical attention for eyes. Do not take internally.

**Shading System**

Flow Plus is available in 16 Vita ccc shades and 5 special shades as follows:

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, Extra Light, G2 (gingival), Universal Opaque, White Opaque, Incisal,

**Notes:**

1. Replace syringe cap after each use to prevent the resin from hardening in the syringe.
2. Flow Plus dispensing tips are intended for single patient use only.

**Storage Conditions**

Store FLOW PLUS at room temperature. It has a shelf life of 2 years from date of manufacture.

(Note: VITA® is a registered trademark of VITA Zahnfabrik Company, Germany).

CAUTION: Uncured resin materials may cause contact dermatitis. Avoid prolonged exposure of uncured resins to skin and soft tissue.

**LIMITED WARRANTY**

**LIMITATION OF MEDICEPT'S LIABILITY**

Medicept's technical advice, whether verbal or in writing, is designed to assist dentists in using Medicept's product. The dentist assumes all risk and liability for damages arising out of the improper use of Medicept's product. In the event of a defect in material or workmanship, Medicept's liability is limited, at Medicept's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Medicept. In no event shall Medicept be liable for any indirect, incidental, or consequential damages.

**EXCEPT AS EXPRESSLY PRO-VIDED ABOVE, THERE ARE NO WARRANTIES, BY MEDICEPT, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.**

**THIS PRODUCT IS SPECIALLY FORMULATED FOR USE IN DENTISTRY**

**Medicept UK Ltd.**

2<sup>nd</sup> Floor, Hygeia House, 66 College Road, Harrow, Middlesex, HA1 1BE, United Kingdom.

Medicept Dental operates a policy of continuous monitoring and improvement of our products. If you have any comments about this product, please contact us at the above address stating the batch number shown on the packaging.